

ZNF134 Antibody



PACO01724

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

ZNF134

Uniprot

P52741

Synonyms:

ZNF134; Zinc finger protein 134


Immunogen:

Synthesized peptide derived from the Internal region of human ZNF134.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using ZNF134 Polyclonal Antibody.