

ACADM Antibody



PACO01743

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ACADM

Uniprot

P11310

Synonyms:

ACADM; Medium-chain specific acyl-CoA dehydrogenase, mitochondrial; MCAD


Immunogen:

Synthesized peptide derived from the Internal region of human MCAD.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using MCAD Polyclonal Antibody.