

VAV1 Antibody



PACO01760

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

VAV1

Uniprot

P15498

Synonyms:

VAV1; VAV; Proto-oncogene vav


Immunogen:

Synthesized peptide derived from human Vav around the non-phosphorylation site of Y160.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using VAV1 Polyclonal Antibody.