

# DFNA5 Antibody



PACO01765

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

GSDME

**Uniprot**

O60443

**Synonyms:**

DFNA5; ICERE1; Non-syndromic hearing impairment protein 5; Inversely correlated with estrogen receptor expression 1; ICERE-1

**Immunogen:**


Synthesized peptide derived from the Internal region of human DFNA5.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of rat stomach 293 cells using DFNA5 Polyclonal Antibody.