

NVL Antibody



PACO01766

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

NVL

Uniprot

O15381

Synonyms:

NVL; Nuclear valosin-containing protein-like; NVLp; Nuclear VCP-like protein


Immunogen:

Synthesized peptide derived from the Internal region of human NVL.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of rat stomach 293 cells using NVL Polyclonal Antibody.