

STAT2 Antibody



PACO01769

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

STAT2

Uniprot

P52630

Synonyms:

STAT2; Signal transducer and activator of transcription 2; p113


Immunogen:

Synthesized peptide derived from human Stat2 around the non-phosphorylation site of Y631.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using Stat2 Polyclonal Antibody.