

TUBA1A/TUBA1B Antibody



PACO01770

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

TUBA1A/TUBA1B

Uniprot

Q71U36/P68363

Synonyms:

TUBA1A; TUBA3; Tubulin alpha-1A chain; Alpha-tubulin 3; Tubulin B-alpha-1; Tubulin alpha-3 chain; TUBA1B; Tubulin alpha-1B chain; Alpha-tubulin ubiquitous; Tubulin K-alpha-1; Tubulin alpha-ubiquitous chain


Immunogen:

Synthesized peptide derived from the C-terminal region of human Tubulin alpha.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% sodium azide.

Product Images



Western Blot analysis of Jurkat k562 cells using Tubulin alpha Polyclonal Antibody.