
Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

GAPDH

Uniprot

P04406

Synonyms:

GAPDH; GAPD; CDABP0047; OK/SW-cl.12; Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH


Immunogen:

Synthesized peptide derived from the N-terminal region of human GAPDH.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.04% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 HepG2 MCF7 cells using GAPDH Polyclonal Antibody.