

SOX17 Antibody



PACO01774

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

SOX17

Uniprot

Q9H6I2

Synonyms:

SOX17; Transcription factor SOX-17


Immunogen:

Synthesized peptide derived from the Internal region of human Sox-17.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Sox-17 Polyclonal Antibody.