

PAOX Antibody



PACO01816

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

PAOX

Uniprot

Q6QHF9

Synonyms:

PAOX; PAO; Peroxisomal N(1)-acetyl-spermine/spermidine oxidase; Polyamine oxidase


Immunogen:

Synthesized peptide derived from the Internal region of human PAOX.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of rat stomach K562 cells using PAOX Polyclonal Antibody.