

ACTC1 Antibody



PACO01828

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ACTC1

Uniprot

P68032

Synonyms:

ACTC1; ACTC; Actin, alpha cardiac muscle 1; Alpha-cardiac actin


Immunogen:

Synthesized peptide derived from the N-terminal region of human Actin-alpha cardiac muscle.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of mouse heart mouse lung cells using Actin-alpha cardiac muscle Polyclonal Antibody.