

TXNDC17 Antibody



PACO01845

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TXNDC17

Uniprot

Q9BRA2

Synonyms:

TXNDC17; TXNL5; Thioredoxin domain-containing protein 17; 14 kDa thioredoxin-related protein; TRP14; Protein 42-9-9; Thioredoxin-like protein 5


Immunogen:

Synthesized peptide derived from the Internal region of human TRP14.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HL-60 Raji K562 SW480 A549 22RV-1 cells using TRP14 Polyclonal Antibody.