

CCL26 Antibody



PACO01864

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CCL26

Uniprot

Q9Y258

Synonyms:

CCL26; SCYA26; C-C motif chemokine 26; CC chemokine IMAC; Eotaxin-3; Macrophage inflammatory protein 4-alpha; MIP-4-alpha; Small-inducible cytokine A26; Thymic stroma chemokine-1; TSC-1


Immunogen:

Synthesized peptide derived from the Internal region of human Eotaxin-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of PC-3 A549 OVCAR3 mouse colon cells using Eotaxin-3 Polyclonal Antibody.