

TFPI2 Antibody



PACO01873

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TFPI2

Uniprot

P48307

Synonyms:

TFPI2; Tissue factor pathway inhibitor 2; TFPI-2; Placental protein 5; PP5


Immunogen:

Synthesized peptide derived from the C-terminal region of human TFPI-2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of H460 293 cells using TFPI-2 Polyclonal Antibody.