

TNFSF13B Antibody



PACO01880

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TNFSF13B

Uniprot

Q9Y275

Synonyms:

TNFSF13B; BAFF; BLYS; TALL1; TNFSF20; ZTNF4; Tumor necrosis factor ligand superfamily member 13B; B lymphocyte stimulator; BLYS; B-cell-activating factor; BAFF;


Immunogen:

Synthesized peptide derived from the Internal region of human TALL-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 22RV-1 HepG2 mouse spleen mouse liver cells using TALL-1 Polyclonal Antibody.