

FGF5 Antibody



PACO01903

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

FGF5

Uniprot

P12034

Synonyms:

FGF5; Fibroblast growth factor 5; FGF-5; Heparin-binding growth factor 5; HBGF-5; Smag-82


Immunogen:

Synthesized peptide derived from the C-terminal region of human FGF-5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of mouse liver cells using FGF-5 Polyclonal Antibody.