

F3 Antibody



PACO01907

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

F3

Uniprot

P13726

Synonyms:

F3; Tissue factor; TF; Coagulation factor III; Thromboplastin; CD142


Immunogen:

Synthesized peptide derived from the Internal region of human TF.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of BXPC-3 mouse eye mouse lung cells using TF Polyclonal Antibody.