

CRP Antibody



PACO01909

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

CRP

Uniprot

P02741

Synonyms:

CRP; PTX1; C-reactive protein


Immunogen:

Synthesized peptide derived from the Internal region of human CRP.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of mouse liver cells using CRP Polyclonal Antibody.