

CDK4 Antibody



PACO01911

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CDK4

Uniprot

P11802

Synonyms:

CDK4; Cyclin-dependent kinase 4; Cell division protein kinase 4; PSK-J3


Immunogen:

Synthesized peptide derived from the N-terminal region of human Cdk4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MCF7 K562 cells using Cdk4 Polyclonal Antibody.