

FSCN1 Antibody



PACO01925

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

FSCN1

Uniprot

Q16658

Synonyms:

FSCN1; FAN1; HSN; SNL; Fascin; 55 kDa actin-bundling protein; Singed-like protein; p55


Immunogen:

Synthesized peptide derived from the Internal region of human Fascin 1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 HuvEc HT29 K562 L929 cells using Fascin 1 Polyclonal Antibody.