

NRG2 Antibody



PACO01927

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

NRG2

Uniprot

O14511

Synonyms:

NRG2; NTAK; Pro-neuregulin-2, membrane-bound isoform; Pro-NRG2


Immunogen:

Synthesized peptide derived from the Internal region of human Neuregulin-2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HT29 cells using Neuregulin-2 Polyclonal Antibody.