

PACO01934

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

LAYN

**Uniprot**

Q6UX15

**Synonyms:**

LAYN; Layilin

**Immunogen:**


Synthesized peptide derived from the Internal region of human Layilin.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of HuvEc K562 cells using Layilin Polyclonal Antibody.