

ECI1 Antibody



PACO01940

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ECI1

Uniprot

P42126

Synonyms:

ECI1; DCI; Enoyl-CoA delta isomerase 1, mitochondrial; 3,2-trans-enoil-CoA isomerase; Delta(3), Delta(2)-enoil-CoA isomerase; D3, D2-enoil-CoA isomerase; Dodecenoyl-CoA isomerase


Immunogen:

Synthesized peptide derived from the C-terminal region of human DCI.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using DCI Polyclonal Antibody.