

ICAM3 Antibody



PACO01956

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

ICAM3

Uniprot

P32942

Synonyms:

ICAM3; Intercellular adhesion molecule 3; ICAM-3; CDw50; ICAM-R; CD50


Immunogen:

Synthesized peptide derived from the Internal region of human CD50.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using CD50 Polyclonal Antibody.