

SLAMF1 Antibody



PACO01980

Product Information

Size:

50ug

Reactivity:

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

SLAMF1

Uniprot

Q13291

Synonyms:

SLAMF1; SLAM; Signaling lymphocytic activation molecule; CDw150; IPO-3; CD150


Immunogen:

Synthesized peptide derived from the region of human SLAM.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using SLAM Polyclonal Antibody.