

MACC1 Antibody



PACO01992

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

MACC1

Uniprot

Q6ZN28

Synonyms:

MACC1; Metastasis-associated in colon cancer protein 1; SH3 domain-containing protein 7a5


Immunogen:

Synthesized peptide derived from the Internal region of human MACC1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of L929 cells using MACC1 Polyclonal Antibody.