

ZKSCAN4 Antibody



PACO02013

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ZKSCAN4

Uniprot

Q969J2

Synonyms:

ZKSCAN4; ZNF307; ZNF427; Zinc finger protein with KRAB and SCAN domains 4; P373c6.1; Zinc finger protein 307; Zinc finger protein 427


Immunogen:

Synthesized peptide derived from the Internal region of human ZKSCAN4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MCF7 HepG2 293 mouse brain mouse lung cells using ZKSCAN4 Polyclonal Antibody.