

NUMB Antibody



PACO02017

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

NUMB

Uniprot

P49757

Synonyms:

NUMB; Protein numb homolog; h-Numb; Protein S171


Immunogen:

Synthesized peptide derived from the Internal region of human NUMB.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using NUMB Polyclonal Antibody.