

MYB Antibody



PACO02027

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

MYB

Uniprot

P10242

Synonyms:

MYB; Transcriptional activator Myb; Proto-oncogene c-Myb


Immunogen:

Synthesized peptide derived from human c-Myb around the non-acetylation site of K471.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using c-Myb Polyclonal Antibody.