

TP73 Antibody



PACO02032

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TP73

Uniprot

O15350

Synonyms:

TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related protein


Immunogen:

Synthesized peptide derived from human p73 around the non-acetylation site of K327.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using p73 Polyclonal Antibody.