

MPO Antibody



PACO02039

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

MPO

Uniprot

P05164

Synonyms:

MPO; Myeloperoxidase; MPO


Immunogen:

Synthesized peptide derived from the N-terminal region of human MPO.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using MPO Polyclonal Antibody.