

ANGPT1 Antibody



PACO02053

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ANGPT1

Uniprot

Q15389

Synonyms:

ANGPT1; KIAA0003; Angiotensin-1; ANG-1


Immunogen:

Synthesized peptide derived from the N-terminal region of human Ang-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Ang-1 Polyclonal Antibody.