

GJA5 Antibody



PACO02064

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

GJA5

Uniprot

P36382

Synonyms:

GJA5; Gap junction alpha-5 protein; Connexin-40; Cx40


Immunogen:

Synthesized peptide derived from the N-terminal region of human Connexin-40.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using Connexin-40 Polyclonal Antibody.