

APAF1 Antibody



PACO02065

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

APAF1

Uniprot

O14727

Synonyms:

APAF1; KIAA0413; Apoptotic protease-activating factor 1; APAF-1


Immunogen:

Synthesized peptide derived from the Internal region of human Apaf-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using Apaf-1 Polyclonal Antibody.