

TYMP Antibody



PACO02099

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TYMP

Uniprot

P19971

Synonyms:

TYMP; ECGF1; Thymidine phosphorylase; TP; Gliostatin; Platelet-derived endothelial cell growth factor; PD-ECGF; TdRPase


Immunogen:

Synthesized peptide derived from the N-terminal region of human PD-ECGF.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MCF-7, 293 cells using PD-ECGF Polyclonal Antibody.