

# UBA52/RPS27A/UBB/UBC Antibody



PACO02111

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

UBA52/RPS27A/UBB/UBC

**Uniprot**

P62987/P62979/P0CG47/P0CG48

**Synonyms:**

UBB; Polyubiquitin-B; UBC; Polyubiquitin-C; RPS27A; UBA80; UBCEP1; Ubiquitin-40S ribosomal protein S27a; Ubiquitin carboxyl extension protein 80; UBA52; UBCEP2

**Immunogen:**


Synthesized peptide derived from human Ub around the non-acetylation site of K63.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of HepG2 cells using Ub Polyclonal Antibody.