

FOXP3 Antibody



PACO02122

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

FOXP3

Uniprot

Q9BZS1

Synonyms:

FOXP3; IPEX; JM2; Forkhead box protein P3; Scurfin


Immunogen:

Synthesized peptide derived from the C-terminal region of human FOXP3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of PC12 cells using FOXP3 Polyclonal Antibody.