

HLA-DQB1/HLA-DQB2 Antibody



PACO02167

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

HLA-DQB1/HLA-DQB2

Uniprot

P01920/P05538

Synonyms:

HLA-DQB1; HLA-DQB; HLA class II histocompatibility antigen, DQ beta 1 chain; MHC class II antigen DQB1; HLA-DQB2; HLA-DXB; HLA class II histocompatibility antigen, DQ beta 2 chain


Immunogen:

Synthesized peptide derived from the Internal region of human HLA-DQB1/2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of PC12 cells using HLA-DQB1/2 Polyclonal Antibody.