

EEF1A1/EEF1A2/EEF1A1P5 Antibody



PACO02169

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

EEF1A1/EEF1A2/EEF1A1P5

Uniprot

P68104/Q05639/Q5VTE0

Synonyms:

EEF1A1; EEF1A; EF1A; LENG7; Elongation factor 1-alpha 1; EF-1-alpha-1; Elongation factor Tu; EF-Tu; Eukaryotic elongation factor 1 A-1; eEF1A-1; Leukocyte receptor cluster member 7; EEF1A2; EEF1AL; STN; Elongation factor 1-alpha 2; EF-1-alpha-2; Eukaryoti


Immunogen:

Synthesized peptide derived from the N-terminal region of human EF-1 alpha1/2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using EF-1 alpha1/2 Polyclonal Antibody.