

UBA52/RPS27A/UBB/UBC Antibody



PACO02172

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

UBA52/RPS27A/UBB/UBC

Uniprot

P62987/P62979/P0CG47/P0CG48

Synonyms:

UBB; Polyubiquitin-B; UBC; Polyubiquitin-C; RPS27A; UBA80; UBCEP1; Ubiquitin-40S ribosomal protein S27a; Ubiquitin carboxyl extension protein 80; UBA52; UBCEP2; Ubiquitin-60S ribosomal protein L40; CEP52; Ubiquitin A-52 residue ribosomal protein fusion pr

Immunogen:

Synthesized peptide derived from the N-terminal region of human Ub.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HepG2 cells using Ub Polyclonal Antibody.

