

NTN1 Antibody



PACO02190

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

NTN1

Uniprot

O95631

Synonyms:

NTN1; NTN1L; Netrin-1; Epididymis tissue protein Li 131P


Immunogen:

Synthesized peptide derived from the C-terminal region of human Netrin-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Netrin-1 Polyclonal Antibody.