

NTS Antibody



PACO02202

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

NTS

Uniprot

P30990

Synonyms:

NTS; Neurotensin/neuromedin N


Immunogen:

Synthesized peptide derived from the C-terminal region of human Neurotensin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA cells using Neurotensin Polyclonal Antibody.