

PACO02230

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

NGF

Uniprot

P01138

Synonyms:

NGF; NGFB; Beta-nerve growth factor; Beta-NGF


Immunogen:

Synthesized peptide derived from the Internal region of human NGF.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa HepG2 cells using NGF Polyclonal Antibody.