

AICDA Antibody



PACO02231

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

AICDA

Uniprot

Q9GZX7

Synonyms:

AICDA; AID; Activation-induced cytidine deaminase; Cytidine aminohydrolase


Immunogen:

Synthesized peptide derived from the Internal region of human AID.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa HepG2 cells using AID Polyclonal Antibody.