

# FCAR Antibody



PACO02233

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

FCAR

**Uniprot**

P24071

**Synonyms:**

FCAR; CD89; Immunoglobulin alpha Fc receptor; IgA Fc receptor; CD89

**Immunogen:**


Synthesized peptide derived from the Internal region of human CD89.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of NIH-3T3 cells using CD89 Polyclonal Antibody.