

SELPLG Antibody



PACO02237

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

SELPLG

Uniprot

Q14242

Synonyms:

SELPLG; P-selectin glycoprotein ligand 1; PSGL-1; Selectin P ligand; CD162


Immunogen:

Synthesized peptide derived from the N-terminal region of human PSGL-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using PSGL-1 Polyclonal Antibody.