

SOD2 Antibody



PACO02240

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

SOD2

Uniprot

P04179

Synonyms:

SOD2; Superoxide dismutase [Mn], mitochondrial


Immunogen:

Synthesized peptide derived from the Internal region of human SOD-2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using SOD-2 Polyclonal Antibody.