

AREG Antibody



PACO02268

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

AREG

Uniprot

P15514

Synonyms:

AREG; SDGF; AREGB; Amphiregulin; AR; Colorectum cell-derived growth factor; CRDGF


Immunogen:

Synthesized peptide derived from the Internal region of human Amphiregulin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HBE cells using Amphiregulin Polyclonal Antibody.