

Phospho-TNNI3 (S22/S23) Antibody



PACO02275

Product Information

Size:

50ug

Reactivity:

Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

TNNI3

Uniprot

P19429

Synonyms:

TNNI3; TNNC1; Troponin I; cardiac muscle; Cardiac troponin I


Immunogen:

Synthesized peptide derived from human Troponin I-C around the phosphorylation site of S22/S23.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of SH-SY5Y cells using Phospho-Troponin I-C (S22/S23) Polyclonal Antibody.