

LAMC3 Antibody



PACO02311

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

LAMC3

Uniprot

Q9Y6N6

Synonyms:

LAMC3; Laminin subunit gamma-3; Laminin-12 subunit gamma; Laminin-14 subunit gamma; Laminin-15 subunit gamma


Immunogen:

Synthesized peptide derived from the C-terminal region of human Laminin gamma-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of SH-SY5Y cells using Laminin gamma-3 Polyclonal Antibody.